

[Continue](#)

Copyright © 2022 Independent Dental • 501 South Beech St • Sutherland • Iowa • 51058 • (800) 728-0020 Site Developed By: Emagine, LLC SDS Kerr OptiMix 100 Operating Guide The OptiMix amalgamator is a fully computerized and programmable 3-touch system used for the mixing of dental alloys, capsule-delivered temporary filling materials, and oral cements. Keypad entry programming Computerized operating system with digital technology LCD display of programmed settings Pre-programmed settings for standard capsules Capsule holder designed for easy capsule loading and removal Safety switch prevents mixing while mixing chamber lid is raised. Internal circuitry clearly labeled with "quick-connect" electrical connections Multi-lingual LCD readout Produces consistent mixes of materials Easy to set mixing frequency and time Quietest triturator Capsule placement and removal can be accomplished with one hand Compatible with all tested capsule designs Frequencies are accurate No undue vibration during operation Technical Specifications: Power requirements: 115/120V~0.5 A, 50/60 Hz, 230V~0.25A, 50/60 Hz ISO 7468 amalgamator class: Type 2 Variable Frequency Dimensions (inches): 8 x 10 x 5 As-Is. No Returns. Packaging: Please allow 1 business day for careful packaging. Shipping: Will ship worldwide! Payment: Payment required within 3 days. Texas residents will have to pay 8.25% sales tax. *The sale of this item may be subject to strict regulation by the U.S. Food and Drug Administration and state and local regulatory agencies. If so, do not bid on this item unless you are an authorized purchaser. If the item is subject to FDA regulation, I will verify your status as an authorized purchaser of this item before shipping of the item.' DISCLAIMER: Regardless of the origin of the equipment, documentation provided or identification appearing upon the equipment, the equipment described and offered here is in no way certified for, recommended for, or offered for any specific use. The purchaser agrees that the seller shall not be held responsible or liable for any injuries or damages, whether incidental or consequential, associated in any way with the equipment. The purchaser, by bidding on this equipment, indicates their acknowledgment of, and agreement to the terms of this disclaimer. OptiBond™ eXTRa UniversalOptiBond™ UniversalOptiBond™ Solo PlusOptiBond™ FL Shape Fill / Obturation Access Diagnose Clean Ultrasonics Retreat ZenFlex™ TF™ Adaptive elements™ e-motionK3™ Nickel-Titanium Files NTI® DiamondsBluWhite Diamonds™ CaviCide™ & CaviWipes™ CaviCide1™ & CaviWipes1™ CaviWipes™ 2.0CaviWipes™ HP belleWax™ Dipping WaxbelleWax™ Esthetic WaxbelleWax™ Sculptring WaxOcclusal Indicator Wax